

## II. Isolate DNA from Plant or Animal Tissue

### REAGENTS, SUPPLIES, & EQUIPMENT

For each group

Container with cracked or crushed ice

DNA rehydration solution (250  $\mu$ L)

70% Ethanol (1.5 mL)

Isopropanol (1.5 mL)

4 microcentrifuge tubes (1.5 mL)

Micropipettes and tips (100-1000  $\mu$ L)

Nuclei Lysis solution (1.5 mL)\*

Permanent marker

Protein Precipitation solution (0.5 mL)

RNAse solution (10  $\mu$ L)

2 Plastic pestles

2 Specimen tissue samples (from Part I)

To share

Microcentrifuge

Water bath or heating block at 65°C

Vortexer (optional)

Microcentrifuge tube rack

\*Store on ice

This universal DNA extraction method uses a commercial kit. Although it is more expensive than the alternate method for plants using Edward's buffer (see part IIa), it has the advantage of working reproducibly with almost any kind of plant or animal specimen.

The large end of a 1000  $\mu$ L pipette tip will punch leaf disks of this size. Animal tissue should be about the size of a pencil eraser. Using more than the recommended amount can inhibit the DNA extraction or amplification.

Lysis solution dissolves membrane bound organelles including the nucleus, mitochondria, and chloroplast.

Grinding the plant tissue breaks up the cell walls. When fully ground, the sample should be a green liquid. There may be some particulate matter remaining.

Step 7 degrades RNA that could interfere with PCR.

Step 9 causes many proteins to precipitate out of the solution, leaving DNA in the supernatant.

1. Obtain plant or animal tissue ~10-20 mg or ¼ inch diameter from your sample. If you are working with more than one sample, be careful not to cross contaminate specimens. (If you only have one specimen, make a balance tube with the appropriate volume of water for centrifuge steps.)
2. Place sample in a clean 1.5 mL tube labeled with an identification number.
3. Add 100  $\mu$ L of nuclei lysis solution to tube.
4. Twist a clean plastic pestle against the inner surface of 1.5 mL tube to *forcefully* grind the tissue for 1 minute. Use a clean pestle for each tube if you are doing more than one sample.
5. Add 500  $\mu$ L more nuclei lysis solution to tube.
6. Incubate the tube in a water bath or heat block at 65°C for 15 minutes.
7. Add 3  $\mu$ L of RNAse solution to tube. Close cap, and mix by rapidly inverting tube several times.
8. Incubate the tube in a water bath or heat block at 37°C for 15 minutes. Then stand tube at room temperature for 5 minutes.
9. Add 200  $\mu$ L of protein precipitation solution to each tube. Vortex tubes for 5 seconds: Securely grasp the upper part of tube, and vigorously hit the bottom end with the index finger of the opposite hand. Use a vortexer if available.
10. Stand tube on ice for 5 minutes.

Centrifugation pellets the nucleic acids. The pellet may appear as a tiny teardrop-shaped smear or particles on the bottom side of the tube underneath the hinge. Do not be concerned if you cannot see a pellet. A large or greenish pellet results when cellular debris carried over from the first centrifugation.

Dry the pellets quickly with a hair dryer. To prevent blowing the pellet away, direct the air across the tube mouth, not into the tube.

In Part III, you will use 2.5  $\mu\text{L}$  of DNA for each PCR. This is a crude DNA extract and contains nucleases that will eventually fragment the DNA at room temperature. Keep the sample cold to limit this activity.

11. Place your tube and those of other groups in a balanced configuration in a microcentrifuge, with cap hinges pointing outward. Centrifuge for 4 minutes at maximum speed to pellet protein and cell debris.
12. Label a clean 1.5 mL tube with your sample number. Use a fresh tip to transfer 600  $\mu\text{L}$  of supernatant to the clean tube. Be careful not to disturb the pelleted debris when transferring the supernatant. Discard old tube containing the precipitate.
13. Add 600  $\mu\text{L}$  of isopropanol to the supernatant in tube. Close cap, and mix by rapidly inverting tubes several times.
14. Place your tube and those of other groups in a balanced configuration in a microcentrifuge, with cap hinges pointing outward. Centrifuge for 1 minute at maximum speed to pellet the DNA.
15. Carefully pour off the supernatant from tube, and add 600  $\mu\text{L}$  of 70% ethanol. Close cap, and flick the bottom of each tube several times to "wash" the pellet.
16. Centrifuge the tube for 1 minute at maximum speed.
17. Carefully pour off the supernatant. Use a micropipette with fresh tip to remove any remaining ethanol, being careful not to disturb the pellet.
18. Air dry the pellet for 10-15 minutes to evaporate remaining ethanol.
19. Add 100  $\mu\text{L}$  of the DNA rehydration solution to each tube, and dissolve the DNA pellet by pipetting in and out several times.
20. Incubate the DNA at 65°C for 45-60 minutes, or overnight at 4°C.
21. Store your sample on ice or at -20°C until you are ready to begin Part III.

## Ila. Isolate DNA from Plant Tissue (Alternate)

### REAGENTS, SUPPLIES, & EQUIPMENT

For each group

Container with cracked or crushed ice

Edward's buffer (2.2 mL)

70% Ethanol

Isopropanol

4 microcentrifuge tubes (1.5 mL)

Micropipettes and tips (100-1000  $\mu\text{L}$ )

Permanent marker

Plant specimens

2 Plastic pestles

3M Sodium Acetate (150  $\mu\text{L}$ )

Tris/EDTA (TE) buffer with RNase A (250  $\mu\text{L}$ )

dH<sub>2</sub>O (1.5 mL)

To share

Microcentrifuge

Vortexer (optional)

Water bath or heating block

This method is optimized for plants. Although it takes about 20 minutes longer than the previous method, it uses readily available reagents.

1. Obtain plant tissue  $\sim\frac{1}{4}$  inch diameter from a specimen. If you are working with more than one sample, be careful not to cross-contaminate the specimens. If you only have one specimen, make a duplicate prep to provide a balance for centrifuge steps.

The large end of a 1000- $\mu\text{L}$  pipette tip will punch disks of this size.

Centrifugation pellets the nucleic acids. The pellet may appear as a tiny teardrop-shaped smear or particles on the bottom side of the tube underneath the hinge. Do not be concerned if you cannot see a pellet. A large or greenish pellet results when cellular debris carried over from the first centrifugation.

Dry the pellets quickly with a hair dryer. To prevent blowing the pellet away, direct the air across the tube mouth, not into the tube.

In Part III, you will use 2.5  $\mu$ L of DNA for each PCR. This is a crude DNA extract and contains nucleases that will eventually fragment the DNA at room temperature. Keep the sample cold to limit this activity.

11. Place your tube and those of other groups in a balanced configuration in a microcentrifuge, with cap hinges pointing outward. Centrifuge for 4 minutes at maximum speed to pellet protein and cell debris.
12. Label a clean 1.5 mL tube with your sample number. Use a fresh tip to transfer 600  $\mu$ L of supernatant to the clean tube. Be careful not to disturb the pelleted debris when transferring the supernatant. Discard old tube containing the precipitate.
13. Add 600  $\mu$ L of isopropanol to the supernatant in tube. Close cap, and mix by rapidly inverting tubes several times.
14. Place your tube and those of other groups in a balanced configuration in a microcentrifuge, with cap hinges pointing outward. Centrifuge for 1 minute at maximum speed to pellet the DNA.
15. Carefully pour off the supernatant from tube, and add 600  $\mu$ L of 70% ethanol. Close cap, and flick the bottom of each tube several times to "wash" the pellet.
16. Centrifuge the tube for 1 minute at maximum speed.
17. Carefully pour off the supernatant. Use a micropipette with fresh tip to remove any remaining ethanol, being careful not to disturb the pellet.
18. Air dry the pellet for 10-15 minutes to evaporate remaining ethanol.
19. Add 100  $\mu$ L of the DNA rehydration solution to each tube, and dissolve the DNA pellet by pipetting in and out several times.
20. Incubate the DNA at 65°C for 45-60 minutes, or overnight at 4°C.
21. Store your sample on ice or at -20°C until you are ready to begin Part III.

## Ila. Isolate DNA from Plant Tissue (Alternate)

### REAGENTS, SUPPLIES, & EQUIPMENT

For each group

Container with cracked or crushed ice

Edward's buffer (2.2 mL)

70% Ethanol

Isopropanol

4 microcentrifuge tubes (1.5 mL)

Micropipettes and tips (100-1000  $\mu$ L)

Permanent marker

Plant specimens

2 Plastic pestles

3M Sodium Acetate (150  $\mu$ L)

Tris/EDTA (TE) buffer with RNase A (250  $\mu$ L)

dH<sub>2</sub>O (1.5 mL)

To share

Microcentrifuge

Vortexer (optional)

Water bath or heating block

This method is optimized for plants. Although it takes about 20 minutes longer than the previous method, it uses readily available reagents.

The large end of a 1000- $\mu$ L pipette tip will punch disks of this size.

1. Obtain plant tissue  $\sim$ 1/4 inch diameter from a specimen. If you are working with more than one sample, be careful not to cross-contaminate the specimens. If you only have one specimen, make a duplicate prep to provide a balance for centrifuge steps.

Detergent in the Edward's buffer, sodium dodecyl sulfate (SDS), dissolves lipids of the cell membranes.

Step 7 denatures proteins, including enzymes that digest DNA.

Step 8 pellets insoluble material at the bottom of the tube.

Step 9 precipitates nucleic acids, including DNA.

Centrifugation pellets the nucleic acids. The pellet may appear as a tiny teardrop-shaped smear or particles on the bottom side of the tube underneath the hinge. Do not be concerned if you can't see a pellet. A large or greenish pellet results when cellular debris carried over from the first centrifugation.)

Nucleic acid pellets are not soluble in ethanol and will not dissolve during washing.

Dry the pellets quickly with a hair dryer. To prevent blowing the pellet away, direct the air across the tube mouth, not into the tube.

If needed, you may store DNA in TE/RNase solution at  $-20^{\circ}\text{C}$  until ready to continue.

2. Place sample in a clean 1.5 mL tube labeled with an identification number.
3. Add 100  $\mu\text{L}$  of Edward's buffer to tube.
4. Grind the tissue for 1 minute by forcefully twisting a clean plastic pestle against the inner surface of the 1.5 mL tube. If you are doing more than one sample, use a clean pestle for each sample.
5. Add 900  $\mu\text{L}$  more Edward's buffer to tube, and grind briefly to remove tissue from the pestle.
6. Vortex tube for 5 seconds, by hand or machine (if available).
7. Boil sample at  $100^{\circ}\text{C}$  for 5 minutes in a water bath or heating block.
8. Place tube, along with those from other groups, in a balanced configuration in a microcentrifuge, and centrifuge for 2 minutes to pellet any remaining cell debris. Centrifuge longer if there is still unpelleted debris.
9. Label a clean 1.5 mL tube with your sample number. If doing more than one sample, use fresh tips to transfer 350  $\mu\text{L}$  of supernatant for each sample to the appropriate fresh tubes. Be careful not to disturb the pelleted debris when transferring the supernatant. Discard old tube containing the precipitate.
10. Add 400  $\mu\text{L}$  of isopropanol to the supernatant in tube. Close cap and mix by rapidly inverting several times.
11. Stand tube at room temperature for 3 minutes.
12. Place your tube and those of other groups in a balanced configuration in a microcentrifuge, with cap hinges pointing outward. Centrifuge for 5 minutes at maximum speed to pellet the DNA.
13. Carefully pour off the supernatant from tube, and add 500  $\mu\text{L}$  of 70% ethanol. Close cap, and flick the bottom of tube several times to "wash" the pellet.
14. Place the tube in a balanced configuration in a microcentrifuge, and spin for 1 minute. Align tubes in the rotor with the cap hinges pointing outward.
15. Carefully pour off the supernatant from tube. Centrifuge the tube again for 30 seconds to force any remaining ethanol to the bottom.
16. Use a micropipette to carefully remove the remaining ethanol from tube. Be careful not to disturb the pellet.
17. Air dry the pellet for 10 minutes to evaporate remaining ethanol.
18. Add 100  $\mu\text{L}$  of TE/RNaseA buffer to tube. Dissolve the nucleic acid pellet by pipetting in and out. Take care to wash down the side of the tube underneath the hinge, where the pellet formed during centrifugation. Use a fresh tip for each tube if you are doing more than one sample.
19. Incubate TE/RNaseA solution at room temperature for 5 minutes.
20. Add 400  $\mu\text{L}$  of  $\text{dH}_2\text{O}$  to tube.
21. Add 50  $\mu\text{L}$  of sodium acetate to tube. Close cap, and mix by rapidly inverting tube several times.
22. Add 550  $\mu\text{L}$  of isopropanol to tube to precipitate the DNA. Close cap, and mix by rapidly inverting tube several times.

*Nucleic acid pellets are not soluble in ethanol and will not dissolve during washing.*

*Dry pellet quickly with a hair dryer. To prevent blowing the pellet away, direct the air across the tube mouth, not into the tube.*

*In Part III, you will use 2.5  $\mu$ L of DNA for each PCR. This is a crude DNA extract and contains nucleases that will eventually fragment the DNA at room temperature. Keep the sample cold to limit this activity.*

23. Stand tube at room temperature for 3 minutes.
24. Place the tube in a balanced configuration in a microcentrifuge, and spin for 5 minutes. Align tubes in the rotor with the cap hinges pointing outward.
25. Carefully pour off the supernatant from tube, and add 500  $\mu$ L of 70% ethanol. Close cap, and flick the bottom of tube several times to “wash” the pellet.
26. Place the tube in a balanced configuration in a microcentrifuge, and spin for 1 minute. Align tubes in the rotor with the cap hinges pointing outward.
27. Carefully pour off the supernatant from tube. Centrifuge the tube again for 30 seconds to force any remaining ethanol to the bottom.
28. Use a micropipette to carefully remove the remaining ethanol from tube. Be careful not to disturb the pellet.
29. Air dry the pellet for 10 minutes to evaporate remaining ethanol.
30. Add 100 $\mu$ L of dH<sub>2</sub>O to tube, and dissolve the DNA pellet by pipetting in and out several times.
31. Store your sample(s) on ice or at -20°C until you are ready to begin Part III.